



LIQUIBAND[®] Rapid[™]



THERE'S FAST THEN THERE'S RAPID

A new generation of *accelerated* topical skin adhesive with reduced and consistent dry time.



CONTINUED INNOVATION

Controlled application and reduced dry time makes LiquiBand® Rapid™ a better solution for effective wound closure.

Accelerated drying
technology

Easy
to use

High volume fill



**THERE'S FAST
THEN THERE'S RAPID**

OCTYL CYANOACRYLATE ACCELERATED ADHESIVE

- ✔ Strong 2-Octyl cyanoacrylate adhesive with added accelerant for faster dry times



WINGED APPLICATOR

- ✔ Signature Liquiband® winged applicator
- ✔ Applicator wings provide ease for use
- ✔ Safe activation and control throughout application
- ✔ Less force required to actuate than market leading competitors
- ✔ Required no assembly before use

WINGED APPLICATOR FEATURE:

✔ LiquiBand® Rapid™ ✘ Dermabond Advanced® ✘ Exofin®



POROUS FOAM TIP

- ✔ New to market porous foam tip to control the application of the adhesive without dripping
- ✔ Broad or thin lines can be applied due to the intuitive highlighter shape
- ✔ Accelerant is incorporated directly into the tip instead of a separate filter

FOAM TIP FEATURE:

✔ LiquiBand® Rapid™ ✘ Dermabond Advanced® ✘ Exofin®

*All sub products

Key benefits that help improve user experience no matter the setting



CONSISTENT DRYTIME

LiquiBand® Rapid™ offers consistent dry times regardless of application technique due to the precise combination of our chemical initiator and tip technology.



CONTROLLED APPLICATION

LiquiBand® Rapid™ expresses a high viscosity adhesive through an innovative tip that provides drip-free priming and even, consistent application.



FLEXIBLE CLOSURE

LiquiBand® Rapid™ is formulated from 2-Octyl cyanoacrylate adhesive to maximize durability and flexibility during topical skin closure that is comparable to market leading competitors.



EFFECTIVE MICROBIAL BARRIER

LIQUIBAND® RAPID™ IS PROVEN EFFECTIVE AGAINST:

- ✓ E. coli
- ✓ E. cloacae
- ✓ E. aeruginosa
- ✓ S. epidermidis
- ✓ S. aureus (MRSA)
- ✓ S. aureus
- ✓ C. albicans
- ✓ Asp. brasiliensis

LiquiBand® Rapid™ offers protection against gram-positive, gram-negative and fungal microbes. The water-resistant microbial barrier offers protection for as long as the adhesive film remains intact.

APPLICATION GUIDE

Please refer to included IFU for comprehensive instructions and warnings on how to apply adhesive



- Remove the LiquiBand® Rapid™ applicator from the package using sterile technique.
- Holding the applicator in one hand, gently squeeze the wings to apply sufficient pressure to crush the inner glass ampoule (a cracking sound will be heard).



- Point the tip downwards, if not already.
- Squeeze the wings and hold until violet adhesive can be seen moistening the foam tip. Once adhesive is on the tip, it is ready to use.
- Use LiquiBand® Rapid™ immediately after crushing as adhesive won't flow freely after a few minutes.



- Hold the wound edges tight together with gloved fingers or forceps.
- Apply topical skin adhesive sparingly as a thin film along the edges of the wound. Avoid heavy application.
- Apply topically do not allow applicator tip or adhesive into the wound.

Accelerated
octyl adhesive

Consistent
application

High strength
and flexibility

Glossy
appearance



Reference	Size	Packaging
72014030*	0.8g	6 pieces / box

*Reference available in our stock at REY Médical SA

IRM 06 0218

Advanced Medical Solutions Limited, Premier Park, 33, Road One, Winsford Industrial Estate, Winsford, Cheshire, CW7 3RT UK © 2020. LiquiBand®, LiquiBand® Rapid™ and the logos contained herein are trade marks of Advanced Medical Solutions Limited.

YOUR PARTNER IN SURGERY SINCE 2003

INNOVATION AND QUALITY AT THE SERVICE OF SURGERY

+41 (0)22 743 24 10 • contact@reymedical.ch • www.reymedical.ch

The best technology available on the market, at the best price, at the right time