



O.R. POSITIONING STRAPS



- Soft and strong fabric keeps patient comfortable
- Fits all size patients
- Non-conductive and incinerable since straps contain no heavy metals
- Securely and conveniently attached to OR table with retaining ring and velcro pads
- Available in reusable and disposable

Code	Description	Length	Packaging
30210	Arm, disposable	3.81 cm x 81.28 cm	50 pieces / box
30220	Knee & Body, disposable	7.62 cm x 152.4 cm	16 pair (2 pieces / bag)


 ISO 13485 CERTIFIED
 FDA REGISTERED
 

 LITBRST001 R2 2/22/18

REY MEDICAL SA: YOUR PARTNER IN SURGERY

• **QUALITY** • **SECURITY** • **ECONOMY**

Visit our website: www.reymedical.ch or call us: 022 743 24 10