



TAILORED SOLUTIONS FOR HERNIA REPAIR
SOLUTIONS POUR LA HERNIE





Inguinal
Ventral
Umbilical

INGUINAL

Open surgery / Chirurgie ouverte

LICHTENSTEIN




CONCEPT 4DDOME®



Laparoscopy (ie)

TAPP/TEP



Biomesh P1



UMBILICAL

Open surgery / Chirurgie ouverte

EXTRAPERITONEAL INTRAPERITONEAL



VENTRAL


Open surgery /
Chirurgie ouverte

Laparoscopy (ie)

EXTRAPERITONEAL INTRAPERITONEAL



Premium®



INGUINAL

Open surgery / Chirurgie ouverte

LICHTENSTEIN

		Ref.	Size / Taille (cm)
4DMesh® - Lightweight PP + PLLA - Weight after resorption 30 (± 10) g/m²			
	Flap / Rabat	FBIO SR F128	12 x 8
	Precut / Prédécoupé	4DMESH M613	6 x 13,5
		4DMESH F812	8,5 x 12,5
		4DMESH MBES	12 x 5
		4DMESH MBEL*	7 x 13
Premium® -Lightweight PP - Weight 37,8 (± 15) g/m²			
		VCB PM P613U	6 x 13,5
		VCB PM 1515U	15 x 15
		VCB PM 3030U	30 x 30
		VCB PM 1530U*	15 x 30

CONCEPT 4DOME®


4DDOME® - Lightweight PP + PLLA - Weight after resorption 30 (± 5) g/m²		
	4D DOME 24 SR	4DDOME® Ø 2,4 (small) + mesh 9 x 5,2
	4D DOME 30 SR	4DDOME® Ø 3 (medium) + mesh 9 x 5,2
	4D DOME 38 SR	4DDOME® Ø 3,8 (large) + mesh 9 x 5,2

*References available on demand

INGUINAL

Laparoscopy (ie)

TAPP/TEP

4DMesh® - Lightweight PP + PLLA - Weight after resorption 30 (± 10) g/m²			
	Flap / Rabat	4DMESH RABA	11 x 14
	Precut / Prédécoupé	4DMESH 1215	12 x 15
		4DMESH 1317	13 x 17
		4DMESH 1216*	12 x 16,5
		4DMESH 1717*	17 x 17
	Anatomique/ Anatomical	4DMESH PRSR	right / droit 10,5 x 14
		4DMESH PRSL	left / gauche 10,5 x 14
		4DMESH PRLR	right / droit 12 x 15
		4DMESH PRLl	left / gauche 12 x 15
		4DMESH PRXR	right / droit 12 x 17
		4DMESH PRXL	left / gauche 12 x 17

Premium® - Lightweight PP - Weight 37,8 (± 15) g/m²			
	Precut / Prédécoupé	VCB PM 1215U*	12 x 15
		VCB PM 1216U	12 x 16,5

Biomesh P1 - PP -Weight 100 (± 15) g/m²		
New  Left  Right	FBIO P1 3D SR	10,5 x 14
	FBIO P1 3D SL	10,5 x 14
	FBIO P1 3D MR	12 x 15
	FBIO P1 3D ML	12 x 15
	FBIO P1 3D LR	12 x 17
	FBIO P1 3D LL	12 x 17

*References available on demand

UMBILICAL

Open surgery / Chirurgie ouverte

EXTRAPERITONEAL

4DVentral® - Midweight PP + PLLA - Weight after resorption 65 (± 5) g/m²			
	New	4DVENT05R0	Ø 5
	New	4DVENT07R0	Ø 7
	New	4DVENT09R0	Ø 9
		4DVENT12R0	Ø 12

INTRAPERITONEAL

CA.B.S.' air® composite - PP + ePTFE - Weight 326.9 (± 2) g/m²		
	CABSTIR 05 T	Ø 5
	CABSAIR C7	Ø 6,8
	CABSAIR C09	Ø 9

VENTRAL

Open surgery / Chirurgie ouverte



EXTRAPERITONEAL

4DVentral® - Midweight PP + PLLA - Weight after resorption 65 (± 5) g/m²		
	4DVENT05R0	Ø 5
	4DVENT07R0	Ø 7
	4DVENT09R0	Ø 9
	4DVENT12R0	Ø 12
	4DVENT1515	15 x 15
	4DVENT1530*	15 x 30
	4DVENT2025	20 x 25
	4DVENT3030	30 x 30

Premium® - Lightweight PP - Weight 37,8 (± 15) g/m²		
	VCB PM 1515U	15 x 15
	VCB PM 3030U	30 x 30
	VCB PM 1530U*	15 x 30

*References available on demand

INTRAPERITONEAL

Intramesh T1 - PP + ePTFE - Weight 120 (± 10) g/m ²		
	FBIO T1 1515	15 x 15
	FBIO T1 2020*	20 x 20
	FBIO T1 3030	30 x 30
	FBIO T1 3050*	30 x 50
	FBIO T1 1015	10 x 15
	FBIO T1 1520	15 x 20
	FBIO T1 2025	20 x 25
	FBIO T1 R12 X	Ø 12

*References available on demand

The implants from this product range are medical devices manufactured by COUSIN BIOTECH S.A.S. The CE conformity has been carried out by the notified body SGS Belgium NV (CEI639). The management system of COUSIN BIOTECH S.A.S is certified for compliance with ISO 13485 standard. Please read carefully the instructions for use before using the device.

Les produits de cette gamme sont des dispositifs médicaux fabriqués par COUSIN BIOTECH S.A.S. L'évaluation de sa conformité CE a été réalisée par l'organisme notifié SGS Belgium NV (CEI639). Le système qualité de COUSIN BIOTECH S.A.S a été certifié conforme aux exigences de la norme ISO 13485. Lire attentivement la notice d'instructions avant utilisation du produit.

Ref : PLVREFFG09 - Version : 22/04/21

Non-contractual photos and texts. Specifications liable to change without notice. Cousin Biotech S.A.S.
Capital of 340 656 € - 398 460 261 RCS Lille - N°TVA FR 34 398 460 261
Photos et textes non contractuels. Spécifications susceptibles de modification sans préavis. Cousin
Biotech S.A.S. au capital de 340 656 € - 398 460 261 RCS Lille - N°TVA FR 34 398 460 261

IMPLANTS	CLASSIFICATION
CA.B.S.'air composite	II b
CA.B.S.'air light	II b
Intramesh T1	II b
Premium	II b
Biomesch A2	II b
Biomesch P1	II b
4DDOME	III
4DMesh	III
4DVentral	III
CA.B.S.'air semi-resorbable	III

YOUR PARTNER IN SURGERY SINCE 2003

INNOVATION AND QUALITY AT THE SERVICE OF SURGERY

+41 (0)22 743 24 10 • contact@reymedical.ch • www.reymedical.ch

The best technology available on the market, at the best price, at the right time