


**MONOPOLAR LAPAROSCOPIC SCISSORS
WITH ERGONOMIC HANDLE**

COMFORT LINE



Code	Type	Usable Length	Cuting Length	Insulating sleeve length
C1223M352G	Curve (Metzenbaum) coated jaw	33 ± 0.5 cm	14 mm	32.8 cm



Adaptability

The **rubber inserts** of the handle are **re-movable** and allow the surgeon to **adapt the instrument** to the morphology of his hand. This provides **maximum comfort** during the procedure.

The ergonomic design offers great comfort and **high handling reliability**.

Safety

Developed with **high quality components** and a **reliable process**, the design and ergonomics of the Monopolar Laparoscopic Scissors provide **safety** and **ease of use**.

Perfect cutting quality, even after prolonged use of the monopolar function.

Thin, sharply curved jaws ensure **high cutting precision**.


Ready to use

Entirely **single-use** device allowing to avoid constraints linked to sterilization, cleaning, stock management, process,...

Delivered in **double sterilization packaging** with **traceability label**.

TECHNICAL DATA

- Handle: Ergonomic double fingers ABS
- Shaft insulation: Sleeve PTFE
- External diameter: 5 mm
- Rotating system: Rotatable at 360°
- Sterile: yes - Ethylene Oxide
- Packaging: Double PE pouch + Tyvek peelable



YOUR PARTNER IN SURGERY SINCE 2003

INNOVATION AND QUALITY AT THE SERVICE OF SURGERY

+41 (0)22 743 24 10

contact@reymedical.ch

www.reymedical.ch

The best technology available on the market, at the best price, at the right time