

 **LAGIS[®]**

SINGLE USE ENDOSCOPIC INSTRUMENT



ENDOSCOPIC INSTRUMENT
SafeGRAB™



PRODUCT CODE	WORKING LENGTH	DESCRIPTION	Quantity: 6 Each/Box
--------------	----------------	-------------	----------------------



DI-5716-33D	33cm	• Fenestrated Grasping Forceps, double action, fine serrated, 360° rotation dial, with on/off ratchet
DI-5716-45D	45cm	



DI-5731-33D	33cm	• Atraumatic Flat Grasping Forceps, double action, 360° rotation dial, with on/off ratchet
DI-5731-45D	45cm	



DI-5670-33D	33cm	• Endo-Clinch Grasping Forceps, double action, fenestrated, 360° rotation dial, with on/off ratchet
DI-5670-45D	45cm	



DI-5901-33D	33cm	• Babcock Grasping Forceps, double action, 360° rotation dial, with on/off ratchet
DI-5901-45D	45cm	



DI-5719-33D	33cm	• CROCE-OLMI Grasping Forceps, single action, atraumatic, fenestrated, curved, 360° rotation dial, with on/off ratchet
-------------	------	--



Specifications subject to change without notice.

© 2013 LAGIS-L04-C-2018.04 Not made with natural rubber latex.

References available on demand

YOUR PARTNER IN SURGERY SINCE 2003

INNOVATION AND QUALITY AT THE SERVICE OF SURGERY

+41 (0)22 743 24 10 • contact@reymedical.ch • www.reymedical.ch

The best technology available on the market, at the best price, at the right time