





GASTRIC CALIBRATION TUBE

The optimal Bariatric
guiding and testing tool



Easy and safe BM test

Four oblong openings arranged in a staggered pattern are located distal to the tube to facilitate pressurisation of the neogastrium.



Precision

Markings to help the anaesthesiologist to adjust the correct depth and placement of the tube.




With or Without Balloon

Non collapsing smooth silicone tube available with or without calibration balloon.

FEATURES

- **Medical Silicone** grade
- **Ideal rigidity** for optimal guiding
- **Optimal positioning** of aspiration and suction holes
- **Maximal performance** for leakage test
- **Markings** from 20 cm to 70 cm
- **Totally atraumatic** due to simple design
- **Available with distal balloon** for any upper GI procedure
- **38 FR gold standard size** for one size "fits all" procedures (3FR is 1mm)



CODES

	Diameter	Length	Balloon
GCT90038B*	38FR	90 cm	-
GCT90038C*	38FR	90 cm	X
GCT90036B	36FR	90cm	-

*References available in our stock at REY Médical SA

YOUR PARTNER IN SURGERY SINCE 2003

INNOVATION AND QUALITY AT THE SERVICE OF SURGERY

+41 (0)22 743 24 10

contact@reymedical.ch

www.reymedical.ch

The best technology available on the market, at the best price, at the right time