


INSUFFLATION LINE



ADVANTAGES

- Insufflation lines are designed to **withstand high gas flow rates** to allow **rapid, controlled** and **safe** insufflation
- Tubing is **flexible, non-obstructive, lightweight** and **highly resistant to plication**
- Device packaged in **double peelable sterilisation bag** and **traceability label**




Reliability

Extra-flexible tubing for anti-plication. The polycarbonate filter and membrane with high hydrophobic, antibacterial and antiviral characteristics provides an excellent barrier to patient/insufflator cross contamination.

Filtration efficiency:

Bacteria: 99.999%


Particles and viruses: 99.99%



Ease of use

The ultra-light tubing avoids placing too much weight on the trocar or pneumoperitoneum needle fitting.

Length of the tubing is 3m.



Compatibility

The various fittings ensure compatibility with all the insufflators. The wide range of fittings makes it possible to homogenise the equipment. The secure luer-lock fitting fits on all available trocars.

YOUR PARTNER IN SURGERY SINCE 2003

INNOVATION AND QUALITY AT THE SERVICE OF SURGERY

+41 (0)22 743 24 10 • contact@reymedical.ch • www.reymedical.ch

The best technology available on the market, at the best price, at the right time