



## SINGLE USE SPECIMEN BAG

Morcellation Range



Automatic Range




Semi-Automatic Range



# Single Use Specimen Bag

- Three different ranges of tissue retrieval systems with various sizes to meet different needs
- Ripstop fabric resistant to tearing following inadvertent puncture
- Secure drawstring closure for increased safety
- Automatic large opening mouth to facilitate tissue capture
- Coloured tabs aids bag deployment and orientation



## PRODUCT SPECIFICATIONS

PRODUCT TYPE	PRODUCT CODE	TROCAR ID (mm)	MOUTH DIAMETER (mm)	LENGTH (mm)	VOLUME (ml)	NUMBER Per BOX
SEMI-AUTOMATIC RANGE  -S : Long string	LTX-TRB-E50*	5	38	90	50	10
	LTX-TRB-E100-S*	10	60	100	100	20
	LTX-TRB-E200-S*	10	70	145	200	20 or 50
	LTX-TRB-E500-S*	12	102	175	500	20
	LTX-TRB-E750	12	102	235	750	10
	LTX-TRB-E1000	12	120	240	1,000	10
	LTX-TRB-E1200	12	170	170	1,200	10
	LTX-TRB-E2000-S*	12	170	240	2,000	10
	LTX-TRB-E2500	12	170	300	2,500	10
MORCELLATION RANGE	LTX-TRB-M2000*	12	140	350	2,000	10
	LTX-TRB-M4000*	15	165	365	4,000	10
AUTOMATIC RANGE	LTX-TRB-A150	10	65	105	150	10
	LTX-TRB-A200*	10	65	150	200	10
	LTX-TRB-A500	10	85	210	500	10
	LTX-TRB-A1000	12	110	150	1,000	10
	LTX-TRB-A1500*	12	110	200	1,500	10
	LTX-TRB-A2000*	12	110	250	2,000	10

References in stock at REY Médical SA

LTX2019TRB01EN

LANTEX is registered trademark of Suzhou Lantex Medical Technology Co. Ltd.  0123

**REY MEDICAL SA: YOUR PARTNER IN SURGERY**

• **QUALITY** • **SECURITY** • **ECONOMY**

Visit our website: [www.reymedical.ch](http://www.reymedical.ch) or call us: 022 743 24 10