



LIQUIBAND® Optima

TOPICAL SKIN ADHESIVE

LiquiBand® Optima Topical Skin Adhesive's unique, winged applicator allows for **effortless activation** and **controlled application**.

Winged applicator

- ▶ Enables quick and easy adhesive expression

Precision tip

- ▶ Ideal for areas of dense hair growth*

Microbial barrier
in 10 seconds

Visidome reservoir

- ▶ Provides visualization of remaining adhesive

Take the **sting** out of wound closure



Designed for closure of clean, fresh wounds with easily-apposed, tension-free edges




Contains a sterile blend of 90% butyl and 10% octyl cyanoacrylate that can be stored at ambient temperatures

*When applying LiquiBand® Optima in areas of dense hair growth, hair must be parted before applying adhesive.



**Ideal for use in any setting, including
Emergency Departments, Surgery
Centers, and Physician Offices**

- ▶ Sting-free application with **10-second dry time**
- ▶ Acts as a **water-resistant microbial barrier**, allowing patients to shower with physician approval immediately
- ▶ Offers **excellent cosmesis**
- ▶ Suitable for the **closure of acute wounds** including scalp, chin, facial lacerations, and punctures
- ▶ **Sloughs off naturally** within 5-14 days after application

 Achieve tensile strength equal to a 4-0 suture with **just one layer** of LiquiBand[®] Optima

Fast, Precise, Comfortable.

Fast: close wounds in as little as 10 seconds

Precise: ideal applicator for facial and scalp wounds

Comfortable: no-sting adhesive for minimal patient discomfort

 10 SECOND SET TIME	 MICROBIAL BARRIER	 PRECISION APPLICATOR	 NO STING	 AMBIENT STORAGE	 MADE IN THE UK
--	---	--	--	--	--

Product	Volume Per Device	Qty Per Box	Product Code
Optima	0.5g	10	OPT001
Optima Mini	0.35g	10	OPT002

YOUR PARTNER IN SURGERY SINCE 2003

INNOVATION AND QUALITY AT THE SERVICE OF SURGERY

+41 (0)22 743 24 10 • contact@reymedical.ch • www.reymedical.ch

The best technology available on the market, at the best price, at the right time