

## SMOKE EVACUATION PENCIL ON/OFF



Code	Description	Packaging	Blade Length
PS11-CMP01	Smoke evacuation pencil with holster Telescopic Function with locking mechanism Teflon coated blade electrode Debit Regulator On/ Off 3 m tubing with multi-input female connector 7mm, 10mm or 22mm inner dia. 3 m cable with 3-pin plug	20 pcs/ box	50 mm



### Adaptability

The 3-pin 4mm banana plug for generator connection is suitable for all H.F. generators on the market.

The length of the electrical cord and the suction line in 3m secure the positioning of the equipment around the patient.

The 22mm female flexible connector of the suction hose allows a universal connection.



### Security

The design of the suction tips and the suction line keeps the suction characteristics of the smoke constant.

An ON/OFF regulator is available on the PS11-CMP01.

The constant flow suction combined with the generator provides a high level of safety for the staff present in the room.

The cable integrated in the suction line and its rotating system provide an excellent grip for a safe movement.

**YOUR PARTNER IN SURGERY SINCE 2003**

**INNOVATION AND QUALITY AT THE SERVICE OF SURGERY**

**+41 (0)22 743 24 10 • [contact@reymedical.ch](mailto:contact@reymedical.ch) • [www.reymedical.ch](http://www.reymedical.ch)**

*The best technology available on the market, at the best price, at the right time*