

SMOKE EVACUATION PENCIL



Code	Description	Packaging	Blade Length
PS11-K50GY	Smoke evacuation pencil with holster Telescopic Function with locking mechanism Teflon coated blade electrode 3 m tubing with 22 mm female connector 3 m cable with 3-prong plug	20 pcs / box	50 mm


REY MEDICAL SA: YOUR PARTNER IN SURGERY

• **QUALITY** • **SECURITY** • **ECONOMY**

Visit our website: www.reymedical.ch or call us: 022 743 24 10



SMOKE EVACUATION PENCIL ON/ OFF



Code	Description	Packaging	Blade Length
PS11-CMP01	Smoke evacuation pencil with holster Telescopic Function with locking mechanism Debit Regulator On/ Off 3 m tubing with 22 mm female connector 3 m cable with 3-prong plug	20 pcs/ box	50 mm

REY MEDICAL SA: YOUR PARTNER IN SURGERY

• **QUALITY** • **SECURITY** • **ECONOMY**

Visit our website: www.reymedical.ch or call us: 022 743 24 10