

LAGIS®



Safe and Effective Innovation By
Collaborating With Surgeons

TROCAR SYSTEM


TROCAR SYSTEM TM
FineSEAL

REY MEDICAL SA: YOUR PARTNER IN SURGERY

• **QUALITY** • **SECURITY** • **ECONOMY**

Visit our website: www.reymedical.ch or call us: 022 743 24 10

LAGIS CLEARVIEW™ OPTICAL TROCAR




ClearView™ Optical Trocar

Set include 1 Sleeve & 1 ClearView Optical Obturator

CODE	DESCRIPTION	LENGTH
TR-06V1	Stability sleeve, 5mm diameter	78mm
TR-06VL1	Stability sleeve, 5mm diameter	100mm
TRO-12V1	Stability sleeve, 12mm diameter	108mm

IMPORTANT: Please refer to Instructions for Use for complete indications, contraindications, warnings and precautions.

LAGIS DILATING TIP TROCAR



Bladeless Dilating Tip

- Minimize abdominal wall and vessel trauma

Improving Patient Outcomes



Dilating Tip Trocar

Set include 1 Sleeve & 1 Dilating Obturator

CODE	DESCRIPTION	LENGTH
TR-06D1	Stability sleeve, 5mm diameter	78mm


Dilating Tip Trocar Kit

Set include 2 Sleeves & 1 Dilating Obturator

CODE	DESCRIPTION	LENGTH
TR-06D2	Stability sleeve, 5mm diameter	78mm
TR-06DL2	Stability sleeve, 5mm diameter	100mm

IMPORTANT: Please refer to Instructions for Use for complete indications, contraindications, warnings and precautions.

LAGIS BLUNT TIP TROCAR



Blunt Tip Hasson Trocar


CODE	DESCRIPTION	LENGTH
TRO-12B1-AH	Smooth sleeve + Hasson Cone, 12mm	108mm

IMPORTANT: Please refer to Instructions for Use for complete indications, contraindications, warnings and precautions.

LAGIS UNIVERSAL TROCAR SLEEVE

Durable Universal Seal

- No reducer caps needed




Recessed Stopcock Valve

- Reduce unintentional opening/closing during surgery

Funnel Lead-In

- Enable smooth instrument exchange



True-Seal™ General Converter

- No reducer caps needed

Flexible, Durable Seal

- Facilitate specimen removal and ensure constant insufflation

Color-Coded Trocar Size

- No confusion between 11's & 12's

Stability Thread Design

- Great abdominal wall retention



UNIVERSAL TROCAR SLEEVE

CODE	DESCRIPTION	LENGTH
TR-06	Stability sleeve, 5mm diameter	78mm
TR-06L	Stability sleeve, 5mm diameter	100mm
TRO-12	Stability sleeve, 12mm diameter	108mm

IMPORTANT: Please refer to Instructions for Use for complete indications, contraindications, warnings and precautions.